## **NMR Precision Teslameter PT2026**

- EN Installation and Safety Manual**
- DE Installations- und Sicherheitshandbuch**
- ES Manual de Instalación y Seguridad**
- FR Manuel d'installation et de sécurité**
- IT Manuale di Installazione e Sicurezza**

**Version 1.0 (Revision 4.0)**

**March 2017**

**REVISION HISTORY**

v. 1.0 r. 1.0	November 2016	First release
v. 1.0 r. 2.0	December 2016	Corrections for CE certification
v. 1.0 r. 3.0	February 2017	Translations: German, Spanish, French, Italian
v. 1.0 r. 4.0	March 2017	Corrections to translations Modify the title, to better reflect the purpose Modify how we refer to the User's Manual

## CONTENTS

---

<b>1- General.....</b>	<b>1</b>
1-1 Scope.....	1
1-2 Intended Use.....	1
1-3 Manufacturer Contact Information.....	1
<b>2- Equipment Specifications .....</b>	<b>2</b>
2-1 Dimensions .....	2
2-2 Ratings.....	4
2-3 Back Panel Connections .....	6
<b>3- Unpacking.....</b>	<b>8</b>
<b>4- Installation .....</b>	<b>9</b>
<b>5- Operation .....</b>	<b>11</b>
5-1 Start-Up.....	11
5-2 Shut-Down .....	11
5-3 Reset and Forced Power-Down .....	12
<b>6- Maintenance and Service .....</b>	<b>13</b>
6-1 Calibration .....	13
6-2 Upgrades .....	13
6-3 Troubleshooting .....	13
6-4 Repairs.....	14
6-5 Warranty .....	14
<b>7- Disposal .....</b>	<b>15</b>

# **BASICS**

## **1-General**

### **1-1 SCOPE**

This manual describes the procedures for safe installation, operation, maintenance, service and disposal of the Metrolab Precision Teslameter PT2026. The latest version of this document can be downloaded from our web site (see Section 1-3).

The User's Manual, "NMR Precision Teslameter PT2026 User's Manual," describes the installation and operation of the software, as well as technical details concerning the instrument's computer interface and software driver. This document is provided in electronic form with the instrument, and can also be downloaded from our web site (see Section 1-3).

### **1-2 INTENDED USE**

The Precision Teslameter PT2026 is a magnetometer, used to very precisely measure the flux density ("field strength") of strong magnetic fields.

### **1-3 MANUFACTURER CONTACT INFORMATION**

Metrolab Technology SA  
110 Chemin du Pont-du-Centenaire  
1228 Plan-les-Ouates, Switzerland

E-mail: [contacts@metrolab.com](mailto:contacts@metrolab.com)


Web site: [www.metrolab.com](http://www.metrolab.com)


Telephone: +41 22 884 33 11

# **BASICS**

## **2-Equipment Specifications**

### **2-1 DIMENSIONS**

Main unit	210 x 125 x 324 mm
Multiplexer	250 x 80 x 180 mm
Probe, model 1226	16 x 12 x 231 mm 

Probe, model 1426	<p>Electronics: 16 x 12 x 231 mm</p> <p>Measurement head:</p> <ul style="list-style-type: none"> <li>• p sample: 9.2 x 6.2 x 31.5 mm</li> <li>• D sample: 16.2 x 6.2 x 31.5 mm</li> </ul> <p>Cable, electronics – head: 50 cm, Ø 4.3 mm</p> 
Probe cable length	<p>10 m; custom upon request</p> <p>100 m max total length (incl. multiplexers)</p>

**2-2 RATINGS**


Power	55 VA, 100 – 240 VAC, 50-60 Hz See warnings below.
Overvoltage	Accepts temporary overvoltages occurring on the mains supply. Transient overvoltages up to the levels of overvoltage category II.
Fuse	3.15 A (T), 5x20 mm, 250 V See warning below.
Environment	Indoor use; no air inlet (IP 50)
Operating temperature	10 – 40 °C
Storage / transport temperature	Deuterium probes: 0 – 80 °C Other: -25 – 80 °C
Altitude	≤ 2000 m
Relative humidity	Maximum 80 % for temperatures up to 31 °C, decreasing linearly to 50 % relative humidity at 40 °C
Pollution	Pollution degree 2: only non-conductive pollution occurs, except that occasionally a temporary conductivity caused by condensation is to be expected.
Magnetic environment	< 0.2 T (PT2026 Main Unit)
Electromagnetic environment	Equipment intended to be used in an industrial electromagnetic environment, class A. See warning below.


**CAUTION**

- ⇒ Always use a three-prong IEC detachable power cord that is grounded, rated for 250 V – 10 A, and with adequate temperature ratings. The cord should satisfy IEC 60245-1 for rubber insulated cables or IEC 60227-1 for PVC insulated cables, as well as all applicable standards and legislation for your region.
- ⇒ To avoid electric shock, the power cord protective grounding conductor must be connected to ground.
- ⇒ The disconnecting device is the power plug. Ensure that it can be reached in case of emergency.
- ⇒ For continued fire protection, use only the specified fuse – see above.
- ⇒ The PT2026 Main Unit contains ferromagnetic components and may be subject to violent forces in a strong magnetic field gradient.
- ⇒ This is a class A product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

## 2-3 BACK PANEL CONNECTIONS

1	Power	See Section 2-2 for ratings and warnings
2	Fuse	See Section 2-2 for ratings and warnings
3	Power switch	<p>The front- and back-panel switches have the same functions:</p> <ul style="list-style-type: none"> <li>• Power on: with power off, hold ~1 sec</li> <li>• Power off: with power on, press momentarily</li> <li>• Factory reset: with power on, hold until power-off (~5 sec)</li> <li>• Forced shutdown: when system is hung, hold until power-off (~10 sec)</li> </ul>
4	Probe	The same connector is used for multiplexers and other PT2026 peripherals; proprietary pin-out.
5	Trigger signal	<p>BNC connector, TTL Level; trigger output or input:</p> <ul style="list-style-type: none"> <li>• Trigger output: flags when field attains a given value; B rising, falling or either; pulse or level output; can drive 50Ω.</li> <li>• Trigger input: determines when measurement is taken; rising or falling edge; high impedance input.</li> </ul>
6	Reference clock	<p>BNC connector:</p> <ul style="list-style-type: none"> <li>• Internal reference output: 10 MHz <math>\pm</math> 1 Hz, can drive 50Ω.</li> <li>• External reference input: 10 MHz <math>\pm</math> 30 Hz, 50Ω input. Allowable level +10 to +20 dBm, which corresponds to a 50% duty cycle square wave with peak-peak amplitude between +2.5 and +5 V.</li> </ul>
7	Ethernet	RJ45

8	USB	USB 2.0, Standard B connector
9	Auxiliary 	Mini DIN 5: 1. On/Off (+5 V): short to ground is like pressing power switch; 2. NMR signal output ( $\pm 5$ V); 3. Modulation output ( $\pm 9$ V, continuous-wave only); 4. Sync output (0–5 V, pulsed-wave only); 5. Ground.

**CAUTION**

- ⇒ To avoid potentially dangerous equipment failures, only Metrolab peripherals, designed for the purpose, should be connected to the Probe connector.
- ⇒ To avoid electric shock, only limited-energy circuits (in the sense of EN 61010-1) that conform to the Safety Extra Low Voltage standard (IEC 60364-4-41) should be connected to the Trigger, Clock and Auxiliary connectors.
- ⇒ Treat the probe connectors with care. They have very tight tolerances, and crushing them, or even dropping them on the floor, may deform the shell, making it difficult or impossible to plug in the probe.
- ⇒ When plugging in the probe, make sure the connector is perfectly aligned with the socket; otherwise it will jam.
- ⇒ To unplug the probe, you must grasp the sleeve and pull straight back to release the locking mechanism. Never pull on the cable.
- ⇒ To unplug the female plug of the MUX cable from a MUX6026 multiplexer, press the ring at the base of the socket to release the locking mechanism.

## BASICS

### 3-Unpacking

## BASICS

### 4-Installation

#### CAUTION

- ⇒ Review the warnings in Sections 2-2 and 2-3.
- ⇒ Make sure you have performed the software installation before powering on the PT2026.
- ⇒ If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- ⇒ The safety of any system incorporating the equipment is the responsibility of the assembler of the system.

- Plug the probe and power cables, as well as the USB cable and/or an Ethernet cable, into the back panel:

For details on the back panel, please see Section 2-3.

- If you have a MUX6026 multiplexer, connect it to the PT2026's Probe socket using the MUX cable, and connect the probes to the multiplexer

ports:

- Plug the PT2026 power plug into a wall socket, the USB cable (if applicable) into your computer, and the Ethernet cable (if applicable) into an Ethernet switch connected to your network.
- Install the PT2026 software, as described in the User's Manual.


# **BASICS**

## **5-Operation**

### **5-1 START-UP**

- Power on the PT2026 by pressing the power switch on the front or rear panel (the two are equivalent). If the PT2026 does not power up, try holding the power button a little longer.
- Wait 2.5 minutes (at 25°C ambient) for the internal clock to warm up and stabilize. During the warm-up time, measurements will be out of specification; in addition, the instrument may not lock to an external reference clock, and measurements will be flagged as “questionable.”
- Place the probe in your magnet:

- Start the PT2026 software as described in the User's Manual.

### **5-2 SHUT-DOWN**

- Stop the software as described in the User's Manual.
- Power off the PT2026 by pressing the power switch on the front or rear panel (the two are equivalent).

### 5-3 RESET AND FORCED POWER-DOWN

If things go awry, here are a few things to try. Start with the first – easier – steps, and continue until you run out of things to try. At that point, see Section 6-Maintenance and Service.

- Use the software controls to stop and restart measuring.
- Use the software controls to reconnect the instrument.
- Restart the PT2026 by powering down and back up.
- Use the software controls to restore the factory defaults.
- Keep the power switch depressed while powering down the PT2026. After approximately five seconds, the instrument should restore the factory settings and then power down.
- The power-down reset described above will not work if the instrument firmware has crashed. In this case, continue to hold the power switch, and after approximately ten seconds, the PT2026 will be forced to power down.
- If the PT2026 firmware detects a serious fault during startup, the instrument will reboot in Device Firmware Upgrade (DFU) mode, ready to receive a firmware upgrade. See Section 6-Maintenance and Service.

## ***BASICS***

### ***6-Maintenance and Service***

#### **CAUTION**

⇒ Do not remove covers. Refer servicing to qualified personnel. The contact information for Metrolab can be found in Section 1-3.

#### **6-1 CALIBRATION**

Recalibrating a PT2026 consists of adjusting its internal reference clock, a procedure only to be performed at the Metrolab factory. However, any properly equipped calibration laboratory can certify the calibration by measuring the internal 10 MHz reference clock output on the back panel. In addition, if you use an external 10 MHz reference clock, the PT2026 does not need to be calibrated at all.

All customers are of course free to fix the intervals at which they check the calibration of their instruments, within the context of their quality assurance policy. Metrolab's recommendation is to check the PT2026 calibration once every year.

#### **6-2 UPGRADES**

Metrolab makes improvements and bug fixes for the PT2026 firmware, software and manual available via its website. The Download page of the Metrolab website (see Section 1-3) will always contain the latest versions.

The PT2026 Main Unit is designed such that you can upgrade the firmware in the field. Firmware upgrades are a delicate procedure, as a failure may render the instrument unusable. Metrolab has made every effort to make the process foolproof, but please take your time and follow the instructions provided with the upgrade carefully.

The PT2026 peripherals (probes, multiplexers, etc.) also contain firmware. However, at this time, this firmware cannot be upgraded in the field; peripherals need to be sent back to the factory to be upgraded.

#### **6-3 TROUBLESHOOTING**

Many problems can be diagnosed over the phone or via e-mail. Should you have trouble with your PT2026, you should first contact your local representative.

Please check the distributor list on our web site (see Section 1-3) for the representative closest you.

If your local representative cannot correct the problem, please contact Metrolab (see Section 1-3). Again, we will try to resolve the problem on the spot, and, if that proves to be impossible, arrange for the repair of your instrument.

#### **6-4 REPAIRS**

Please contact Metrolab (see Section 1-3). Repairs should only be performed by the manufacturer. We will issue a Return Material Authorization (RMA) number and provide shipping instructions. If the instrument is out of warranty, we will diagnose the problem and send you an estimate before undertaking the repair.

#### **6-5 WARRANTY**

The standard warranty period is two years from the date of purchase. During this period, Metrolab will repair or replace a failing unit free of charge, unless it is clear that the unit has been abused (crushed probe or electronics, torn cable, etc.). We do not assume responsibility for consequential damage, for example to your PC.

## ***BASICS***

### ***7-Disposal***

#### **NOTICE**

⇒

This product conforms to the WEEE Directive of the European Union (2002/96/EC) and belongs to Category 9 (Monitoring and Control Instruments). For proper environment friendly disposal, you can return the instrument free of charge to us or our local distributor.

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