

Global range of
magnetic measurement
equipments

METROLAB

Instruments SA

- NMR
- ESR
- Hall effect
- Digital fluxmeter
- from mT to 50 T
- Accuracy 10^{-2} to 10^{-6}

tel : +41 (22) 884 33 11 fax : +41 (22) 884 33 10
e-mail: contact@metrolab.com

NEW NMR MAGNETIC
FIELD CAMERA
MFC-3045

for

- MRI magnet mapping
- MRI decay monitoring

- Field range 0.08 T to 7 T
- Mapping accuracy ± 0.2 ppm
- Decay measurement ± 0.03 ppm