




PWP1426 Remote Head for Medium and High Field



PWP1426 Remote Head for D2O sample



Type	L [mm]	W [mm]	T [mm]	P [mm]
Proton sample	31.9 ± 0.5	9.6 ± 0.5	5.9 ± 0.5	4.0 ± 0.3
Deuterium sample	31.9 ± 0.5	16.9 ± 0.5	6.2 ± 0.5	4.0 ± 0.3

A	Initial release	23.05.2025	J. Vieira	J. V.
REV.	DESCRIPTION	DATE	AUTEUR	VERIFIE
Titre:	Remote Head (Proton & Deuterium) NMR Probe 1426		Tolérance gén. $\pm 0.1 \text{ mm}$	Echelle: 2:1
				Format: A4
6		Variant:	No.de dessin:	22255
			Dessin d' ensemble:	Feuille#: 1 / 1